


Enzyme kinetics of alkaline phosphatase



LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

The class

- Alkaline phosphatase kinetics
- Brief background of steady state kinetics and the logic behind the measurements
- Assays – parameters to think about when designing assay, how this assay differs from glucose oxidase assay etc.
- Other assays

LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

Learning outcomes

- What is an enzymatic assay
- Measuring products generated by an enzyme
- Basis of design of assays
- Determination of kinetic parameters

LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

Alkaline phosphatase

- Alkaline phosphatase is an enzyme that hydrolyses esters of phosphoric acid into the corresponding alcohol and phosphoric acid
- Can be isolated from bone, kidney, intestine, plasma, liver, spleen, plants and microorganisms
- The kinetics of enzymes is often studied with model substrates that yield a color, so the progress of the reaction can be followed with a spectrophotometer

LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

Assay nitty-gritties

- How do you quantitate *p*NPP ?

Ans: Calibration curve - graph showing the analytical response as a function of the known quantity of analyte

- What are you calibrating ?

Ans: The color due to reducing sugars

- How ?

Ans: By running standards of reducing sugar (known quantity of analyte) and measuring absorbance (analytical response)

Then plotting the analyte response as a function of reducing sugar

LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

Calibration curve prep


- Here, known is the concentration of sugar you prepare
- Unknown is the absorbance due to enzyme assay
- Record the absorbance of the alcohol
- Plot Abs vs ester
- Use the Abs of *p*NPP generated to determine the concentration or amount of reducing sugar in unknown

LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

Calibration : Best fit

R^2 : compares the sums of the variations for the y-values to the best-fit line relative to the variations to a horizontal line.

- $R^2 \times 100$: percent of the variation of the y-variable that is explained by the variation of the x-variable.
- A perfect fit has an $R^2 = 1$; no relationship for $R^2 \approx 0$




Assay

- First determine the calibration curve by measuring absorbance of the known and plot
- Then measure the concentration of unknown and convert the absorbance to the known to quantify the unknown
- How will you use to determine the MM kinetics ?

The reaction mechanism

LS2102 Biochemistry lab |
Copyright @ Supratim Datta,
IISER Kolkata

Step 1- Note the active site water, hydroxyl ion, Zn and Mg ion.
Hydroxyl ions are formed under alkaline conditions




LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

Comment to previous page:

- In the first step, preceding page. Negatively charged substrate replaces the water and associates with enzyme through the positively charged Arg166 and two Zn^{+2}
- Enzyme substrate complex is formed

LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

Step 2: Deprotonated Ser102 is a strong nucleophile that substitutes *p*-nitrophenol in a nucleophilic substitution reaction. Pi remains covalently bound to Ser102




LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

Comments to prev slide

- You can see that the Zn^{2+} associates with the hydroxyl ion
- The enzyme product complex is formed

LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

Step 3: the hydroxyl ion carries out a nucleophilic substitution reaction at the phospho-serine




LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

Comments to prev slide

- Observe that the bond between the phosphate and Ser102 breaks.
- The inorganic hydroxy phosphate is loosely attached, hence the E dot Pi nomenclature.

LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

Step 4: hydroxy phosphate is replaced by water




LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

Comments to prev slide

- Now the active site is as before and ready for the second reaction. So what is broken step by step is for one turnover (remember k_{cat}). The enzyme catalyzes many such turnover (substrate to product)
- A simple organic chemistry happening at the enzyme active site. Difference is the amino acid side chain spatial arrangements and the cofactors
- Ignoring the formation of ester from alcohol and phosphate is one approximation. If you approximate further by assuming that the ester bond breaking and P_i formation is the rate limiting step, you can simplify the mechanism further to $E+S \rightleftharpoons ES \rightarrow E+P$ MM kinetics

LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata


Steady-state enzyme kinetics



Assumptions of Michaelis-Menten kinetics:

1. The reaction is at equilibrium
2. The reaction is at steady-state

Michaelis-Menten Enzyme Kinetics



Steady-state Model:

- Formation of ES = Disappearance of ES

$$k_1[E][S] = k_2[ES] + k_{\text{cat}}[ES]$$

- Define a new constant, K_m (Michaelis constant):


$$K_m = (k_2 + k_{\text{cat}})/k_1$$

K_m reflects the affinity of the enzyme for the substrate:

- low K_m = high affinity
- $K_m = [S]$ that produces $V = 1/2 V_{\text{max}}$

I S2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

An enzyme's response to substrate can be visualized using the Michaelis-Menten plot



Michaelis-Menten Equation
$$v_o = \frac{V_{\text{max}} [S]}{K_M + [S]}$$

I S2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

Simplifying M & M analysis using a Lineweaver-Burk Plot

- A **double-reciprocal** representation of V vs. [S] data:
plot $1/V$ vs. $1/[S]$
- Take the inverse of both sides of the M-M equation to get the L-B equation, which specifies a **line**

M-M $V = \frac{V_{\max}[S]}{K_m + [S]}$


Inverse M-M $1/V = (K_m + [S]) / V_{\max}[S]$

Rearranging

L-B $1/V = 1/V_{\max} + K_m/V_{\max} (1/[S])$

LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata


Determining Kinetic Parameters from an L-B Plot



$$1/V = 1/V_{\max} + K_m/V_{\max} (1/[S])$$

LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata


Michaelis-Menten Kinetics with Inhibition



Just to give you an idea of how the plot changes in presence of various kinds of inhibitors


LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

Lineweaver-Burk Plots with Inhibitors



LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

Another type of colorimetric reaction



LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata


Why *p*-Nitrophenyl β -D-glucopyranoside and not glucose ?

- Colorimetric assays are simple to perform and easy to measure
- pNP gives a yellow color
- 1 mole of pNP = 1 mole of glucose
- Measure amt of *p*NP => Know amount of glucose
- How will you measure amount of pNP ?

Ans: Calibration curve !

LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

Glucose oxidase test – enzyme based test for glucose





glucose oxidase converts glucose to gluconate and hydrogen peroxide (H_2O_2).

The enzyme peroxidase then catalyzes a redox reaction in the presence of H_2O_2 that changes a colorless reduced dye into a colored oxidized product.

- The level of oxidized dye in the sample is measured using a blood glucose

LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

Quantification of the released phosphate with Malachite Green and ammonium molybdate



LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

Assay

- how much product is formed over a given time or, in some cases, how much substrate has been used up
- a method for measuring either product or substrate in the presence of the other is required.
- **‘stopped assays’** - stopping the reaction after a fixed time, then measuring how much product has been formed.
- chemical, enzymatic to bioassay the simplest is chosen provided it is reliable.
- phosphate release from a phosphate ester can be measured by the standard phosphomolybdate procedure No separation needed.
- Stopping the reaction - denature the enzyme, such as strong acid, alkali or detergent; heat;

LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

When to stop a stopped reaction?


- How do we decide the time of the reaction ?
- these methods provide only one single point and the velocity must be calculated from the slope of a line connecting this point with the blank before starting the reaction.
- control measurements at different reaction times must be undertaken to establish the linearity of the assay and appropriate choice of time point.

LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata

Stopped assay

- Where no direct signal for the conversion of substrate or product can be found, the reaction must be stopped after a defined time and the amount of product formed or substrate converted must be analysed
- Instead of a continuous progress curve these methods provide only one single point and the velocity must be calculated from the slope of a line connecting this point with the blank before starting the reaction.
- control measurements at different reaction times must be undertaken to establish the linearity of the assay and appropriate choice of time point.

LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata



LS2102 Biochemistry lab | Copyright @ Supratim Datta, IISER Kolkata